



Anti-*C. trachomatis* Polyclonal antibody (DPAB-DC4082)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>C. trachomatis</i>
Immunogen	MOMP (major outer membrane protein) from <i>Chlamydia trachomatis</i> .
Source/Host	Goat
Species Reactivity	<i>C. trachomatis</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	1 mg
Buffer	In PBS, pH 7.4 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	<i>Chlamydia trachomatis</i> , an obligate intracellular human pathogen, is one of four bacterial species in the genus <i>Chlamydia</i> . <i>C. trachomatis</i> is a gram-negative bacterium, therefore its cell wall components retain the counter-stain safranin and appear pink und
Keywords	<i>C. trachomatis</i> ; <i>Chlamydia trachomatis</i>