



Anti-Chlormadinone acetate polyclonal antibody (DPAB-DC4084)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Chlormadinone acetate conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Chlormadinone acetate (INN, USAN, BAN, JAN; sold under brand names including Clordion, Gestafortin, Lormin, Non-Ovlon, Normenon, Verton, and many others), sometimes abbreviated as CMA, is a steroidal progestin with additional antiandrogen and antigonadotr
Keywords	CMA; Chlormadinone acetate