



Rabbit Anti-Cholera Toxin Subunit B Polyclonal antibody (DPAB-DC4085)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Beta subunit of Cholera toxin.
Immunogen	Purified cholera toxin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bacteria
Conjugate	Unconjugated
Applications	WB, IHC, IF, IP, ELISA,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Cholera Toxin beta subunit
----------	----------------------------