



Anti-Cholesterol Oxidase (native protein) polyclonal antibody (DPAB-DC4086)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified Microorganism Cholesterol Oxidase.
Source/Host	Goat
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IP, ELISA,
Format	Liquid
Size	1 mg
Buffer	In 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This enzyme belongs to the family of oxidoreductases, specifically those acting on the CH-OH group of donor with oxygen as acceptor. The systematic name of this enzyme class is cholesterol:oxygen oxidoreductase. Other names in common use include cholesterol
Keywords	Cholesterol Oxidase