



Anti-Chromogranin B (C-terminal) polyclonal antibody (DPAB-DC4087)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This sequence is the same in rat, mouse, and human chromogranin B.
Immunogen	C-terminus 15 amino acids of Chromogranin B coupled to ovalbumin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC,
Format	Liquid
Size	200 µl
Buffer	In serum
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Granin (chromogranin and secretogranin) is a protein family of regulated secretory proteins ubiquitously found in the cores of amine and peptide hormone and neurotransmitter dense-core secretory vesicles.
---------------------	--

GENE INFORMATION

Entrez Gene ID [1114](#)

UniProt ID [P05060](#)
