



Anti-ch-TOG (C-terminal) polyclonal antibody (DPAB-DC4088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes human and mouse ch-TOG. A BLAST search demonstrated 100% identity to the peptide immunogen in human, mouse, rat, chimpanzee, dog, cow, chicken, opossum, and frog.
Immunogen	A synthetic peptide corresponding to C-terminus amino acids of ch-TOG.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Chicken, Chimpanzee, Dog, Frog, Human, Mouse, Opossum, Rat
Conjugate	Unconjugated
Applications	WB (Cell lysate), IF, IP,
Format	Liquid
Size	100 µl
Buffer	In antiserum
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords ch-TOG

GENE INFORMATION

Entrez Gene ID [9793](#)

UniProt ID [Q14008](#)
