



# Anti-C. botulinum type A Neurotoxin Polyclonal antibody (DPAB-DC4102)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	C. botulinum Type A Neurotoxin.
Target	C. botulinum type A Neurotoxin
Immunogen	Clostridium botulinum type A Neurotoxin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	C. botulinum
Conjugate	Unconjugated
Applications	WB, ELISA, Dot, EIA,
Format	Liquid
Size	200 µg
Buffer	In PBS, pH 7.2 (20% glycerol)
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Clostridium botulinum is a Gram-positive, rod-shaped, anaerobic, spore-forming, motile
--------------	--

bacterium with the ability to produce the neurotoxin botulinum. The botulinum toxin can cause a severe flaccid paralytic disease in humans and animals and is the most p

---

**Keywords**

Neurotoxin; Clostridium botulinum Type A Neurotoxin

---