



Anti-C. botulinum type B Neurotoxin Polyclonal antibody (DPAB-DC4103)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	C. botulinum Type B Neurotoxin Complex, neurotoxin or protein isolates.
Target	C. botulinum type B Neurotoxin
Immunogen	Clostridium botulinum type B Neurotoxin and complexing proteins.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	C. botulinum
Conjugate	Unconjugated
Applications	WB, ELISA, Dot, EIA,
Format	Liquid
Size	200 µg
Buffer	In PBS, pH 7.2 (20% glycerol)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Clostridium botulinum is a Gram-positive, rod-shaped, anaerobic, spore-forming, motile bacterium with the ability to produce the neurotoxin botulinum. The botulinum toxin can cause a
---------------------	--

severe flaccid paralytic disease in humans and animals and is the most p

Keywords

Neurotoxin Complex; Clostridium botulinum Type B Neurotoxin Complex
