



Anti-Collagen Type I CnBr fragment (native) polyclonal antibody (DPAB-DC4108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It exhibits very moderately reactivity with human and weak with bovine collagen type IX. No cross reactivity with collagen type II.
Immunogen	Native purified collagen Type I CnBr fragment from fetal mouse skin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF, IP, ELISA, EIA, Flow Cyt,
Format	Liquid
Size	100 µg
Buffer	In PBS
Preservative	None
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Collagen Type I CnBr fragment
-----------------	-------------------------------