



Anti-Collagen Type I (native) polyclonal antibody (DPAB-DC4109)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with tuna fish collagen type I.
Immunogen	Native purified Collagen Type I from tuna fish skin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fish
Conjugate	Unconjugated
Applications	IHC-P, IF, ELISA, RIA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from phosphate buffer
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Collagen Type I
-----------------	-----------------