



# Anti-Collagen Type I, II, III, IV, V (native) polyclonal antibody (DPAB-DC4111)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Recognizes human collagen type I, II, III, IV and V.
Immunogen	Native purified human Collagen Type I, II, III, IV, V.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-Fr, ELISA, Dot,
Format	Liquid
Size	100 µl
Buffer	In purified antibody
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Collagen Type I, II, III, IV, V
----------	---------------------------------