



Anti-Collagen Type I/III CnBr fragment (native) polyclonal antibody (DPAB-DC4112)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It exhibits no cross reactivity.
Immunogen	Native purified collagen Type I/III CnBr fragment from fetal mouse skin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF, IP, ELISA, EIA, Flow Cyt,
Format	Liquid
Size	100 µg
Buffer	In PBS
Preservative	None
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Collagen Type I/III CnBr fragment
----------	-----------------------------------