



Rabbit anti-human Collagen Type IV polyclonal antibody (DPAB-DC4120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified human Collagen Type IV.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-Fr,
Format	Liquid
Size	250 µl
Buffer	In Tris-HCl buffer (1% BSA, 0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Collagen IV (ColIV or Col4) is a type of collagen found primarily in the basal lamina. The collagen IV C4 domain at the C-terminus is not removed in post-translational processing, and the fibers link head-to-head, rather than in parallel. Also, collagen I
Keywords	ColIV; Collagen Type IV

GENE INFORMATION

Entrez Gene ID [1282](#)

UniProt ID [P02462](#)
