



Anti-Collagen type V (native) polyclonal antibody (DPAB-DC4121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It exhibits slightly with collagen type IV.
Immunogen	Native purified collagen type V from human placenta.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Marmoset, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF, IP, ELISA, EIA, Flow Cyt,
Format	Liquid
Size	100 µg
Buffer	In PBS
Preservative	None
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Type-V collagen is a form of fibrillar collagen associated with classical Ehlers-Danlos syndrome. Autoimmunity against type V collagen is associated with lung transplant failure.
---------------------	--

Keywords

Collagen type V
