



# Anti-Complement C4c (native) polyclonal antibody (DPAB-DC4125)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native purified complement C4c from human plasma.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Complement C4c
----------	----------------