



# Anti-Corticosterone 21-acetate polyclonal antibody (DPAB-DC4128)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	4-PREGNEN-11b, 21-DIOL-3, 20-DIONE 21-ACETATE 3-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Corticosterone 21-acetate
----------	---------------------------