



Anti-CORT (C-terminal) polyclonal antibody (DPAB-DC4129)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	4-PREGNEN-11beta, 21-DIOL-3, 20-DIONE 3-O-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Corticosterone (CORT) is a 21-carbon steroid hormone of the corticosteroid type produced in the cortex of the adrenal glands in rodents and other non-human animals.
Keywords	CORT; Corticosterone