



Anti-Cortisone (C-terminal) polyclonal antibody (DPAB-DC4133)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Immunogen | 4-PREGNEN-17, 21-DIOL-3, 11, 20-TRIONE 21-ACETATE 3-CARBOXYMETHYLOXIME. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | N/A |
| Conjugate | Unconjugated |
| Applications | IP, EIA, |
| Format | Liquid |
| Size | 150 µg |
| Buffer | In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal) |
| Preservative | None |
| Storage | Store at -20°C. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Cortisone (17-hydroxy-11-dehydrocorticosterone) is a 21-carbon steroid hormone. It is one of the main hormones released by the adrenal gland in response to stress. In chemical structure, it is a corticosteroid closely related to cortisol. It is used to tr |
| Keywords | Cortisone |