



Anti-CTTN polyclonal antibody (DPAB-DC4141)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to human CTTN.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Rat
Purification	Protein A purification
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Concentration	1 mg/mL
Size	100 µl
Buffer	In HEPES, 150 mM NaCl (50% glycerol, 0.01% BSA, 0.03% sodium azide)
Preservative	0.03% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CTTN
-----------------	------

GENE INFORMATION

Entrez Gene ID [2017](#)

UniProt ID [Q14247](#)
