



Anti-CXADR (C-terminal) polyclonal antibody (DPAB-DC4143)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminus of rabbit CXADR.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rabbit
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CXADR
----------	-------

GENE INFORMATION

Entrez Gene ID [1525](#)

UniProt ID [P78310](#)
