



Anti-CXCL10 polyclonal antibody (DPAB-DC4144)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to 83 amino acids of bovine CXCL10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Goat, Sheep
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CXCL10

GENE INFORMATION

Entrez Gene ID [3627](#)

UniProt ID [P02778](#)
