



Anti-Cy3/Cy5 polyclonal antibody (DPAB-DC4145)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reactivity was observed against Cy2, Cy3 and Cy5.
Immunogen	A mixture of Cy3 and Cy5 conjugated with KLH.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC, ELISA,
Format	Liquid
Size	1 mg
Buffer	In 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Cy3/Cy5
----------	---------