



Anti-CYP1A1/CYP1A2 polyclonal antibody (DPAB-DC4147)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This polyclonal goat IgG is useful for identification of CYP450 1A1 and 1A2 forms by immunoassay or Western blot. Reacts with both forms in liver microsomes prepared from rabbits or rats that were induced to express multiple forms of P450s.
Immunogen	Purified rabbit CYP1A1/CYP1A2.
Source/Host	Goat
Species Reactivity	Rabbit, Rat
Conjugate	Unconjugated
Applications	WB, IP,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.05% sodium azide
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CYP1A1/CYP1A2

GENE INFORMATION

Entrez Gene ID [1543](#)

UniProt ID [A0N0X8](#)
