



Anti-CYP2E1 polyclonal antibody (DPAB-DC4148)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This polyclonal goat IgG made against rabbit CYP450 2E1 is useful for identification of this ethanol and pyridine inducible P450 form by immunoassay or Western blots. Reacts as a single band with liver microsomes obtained from induced (untreated and treat
Immunogen	Rabbit CYP2E1.
Source/Host	Goat
Species Reactivity	Rabbit, Rat
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µl
Buffer	In purified antibody
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CYP2E1

GENE INFORMATION

Entrez Gene ID [1571](#)

UniProt ID [P05181](#)
