



Anti-CYP450 2B4 polyclonal antibody (DPAB-DC4149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This goat polyclonal antibody specifically reacts with rabbit CYP450 2B4. Reacts as a single band with liver microsomes prepared from rabbits that were induced to express multiple forms of P450s.
Immunogen	Rabbit CYP450 2B4.
Source/Host	Goat
Species Reactivity	Rabbit, Rat
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.05% sodium azide
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CYP450 2B4