



Anti-Cytomegalovirus AD169 Polyclonal antibody (DPAB-DC4152)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with various antigens from all CMV strains. No cross reactivity with other herpes viruses.
Target	Cytomegalovirus AD169
Immunogen	Cytomegalovirus strain AD169, whole virus.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Cytomegalovirus
Conjugate	Unconjugated
Applications	IF, sELISA, Capture Antibody, Detection Antibody,
Format	Liquid
Size	100 µl
Buffer	In antiserum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Cytomegalovirus is a viral genus of the viral family known as Herpesviridae or herpesviruses. It is typically abbreviated as CMV. The species that infects humans is commonly known as human CMV (HCMV) or human herpesvirus-5 (HHV-5), and is the most studied

Keywords

CMV AD169; Cytomegalovirus (AD 169)
