



Anti-DAP12 polyclonal antibody (DPAB-DC4155)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to DAP12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	IP,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DAP12
----------	-------