



Anti-DEFB103A (N-terminal) polyclonal antibody (DPAB-DC4161)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with beta-defensin-3 of human and mouse origin, N-terminal fragments of beta-defensin-3, no cross-reactivity with other known beta-defensins or alpha-defensins.
Immunogen	A mixture of synthetic peptides (conjugated with KLH) corresponding to C-terminus 15 amino acids and N-terminus 13 human and mouse DEFB103A/MBD-3. The sequence is LPKEEQIGKCSTRGR
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	ELISA, RIA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from purified antibody in Tris buffer containing BSA.
Preservative	None
Storage	Store at -20°C on dry atmosphere. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords DEFB103A

GENE INFORMATION

Entrez Gene ID [414325](#)

UniProt ID [P81534](#)
