



# Sheep anti-Desipramine polyclonal antibody (DPAB-DC4162)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Desipramine conjugated with BTG.
Source/Host	Sheep
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Desipramine (also known as desmethylinipramine) is a tricyclic antidepressant (TCA). It inhibits the reuptake of norepinephrine and to a minor extent serotonin. It is used to treat depression, but not considered a first line treatment since the introduction of selective serotonin reuptake inhibitors (SSRIs).
Keywords	Desipramine