



Anti-Desoxycorticosterone polyclonal antibody (DPAB-DC4163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	4-PREGNEN-21-OL-3, 20-DIONE 3-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	11-Deoxycorticosterone (DOC), or simply deoxycorticosterone, also known as 21-hydroxyprogesterone, is a steroid hormone produced by the adrenal gland that possesses mineralocorticoid activity and acts as a precursor to aldosterone. It has no significant g
Keywords	DOC; Desoxycorticosterone