



# Anti-DEX (native) polyclonal antibody (DPAB-DC4164)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native purified Penicillium DEX .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fungi
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	DEX ; DEX
----------	-----------