



Anti-Dexamethasone polyclonal antibody (DPAB-DC4165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Flumethasone 50%, Betamethasone 43%
Immunogen	Dexamethasone conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Dexamethasone is a potent synthetic member of the glucocorticoid class of steroid drugs that has anti-inflammatory and immunosuppressant effects. It is 25 times more potent than cortisol in its glucocorticoid effect, while having minimal mineralocorticoid
Keywords	Dexamethasone