



Anti-Digoxigenin polyclonal antibody (DPAB-DC4174)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Digoxigenin conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Digoxigenin (DIG) is a steroid found exclusively in the flowers and leaves of the plants <i>Digitalis purpurea</i> , <i>Digitalis orientalis</i> and <i>Digitalis lanata</i> (foxgloves), where it is attached to sugars, to form the glycosides.
Keywords	DIG; Digoxigenin