



Anti-Dihydrotestosterone (C-terminal) polyclonal antibody (DPAB-DC4176)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	5a-ANDROSTAN-17b-OL-3-ONE3-O-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	RIA, EIA
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Dihydrotestosterone (commonly abbreviated to DHT), or 5α-dihydrotestosterone (5α-DHT), also known as androstanolone (5α-androstan-17β-ol-3-one) as well as 17β-hydroxy-5α-androstan-3-one, is a sex steroid and androgen hormone. The enzyme 5α-reductase synth
Keywords	DHT; Dihydrotestosterone