



Anti-Dihydrotestosterone polyclonal antibody (DPAB-DC4178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Dihydrotestosterone conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Dihydrotestosterone (commonly abbreviated to DHT), or 5 α -dihydrotestosterone (5 α -DHT), also known as androstanolone (5 α -androstane-17 β -ol-3-one) as well as 17 β -hydroxy-5 α -androstane-3-one, is a sex steroid and androgen hormone. The enzyme 5 α -reductase synth
Keywords	DHT; Dihydrotestosterone