



Anti-Dimethylbenzanthracene polyclonal antibody (DPAB-DC4179)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DMBA conjugated with BSA.
Source/Host	Rat
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IHC-P, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from antiserum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DMBA; Dimethylbenzanthracene (DMBA)
----------	-------------------------------------