



# Anti-Diphtheria Diphtheria Toxin/Toxoid polyclonal antibody (DPAB-DC4180)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Demonstrated reactivity to purified diphtheria toxin/toxoid.
Immunogen	Purified diphtheria toxoid.
Source/Host	Goat
Species Reactivity	Bacteria
Conjugate	Unconjugated
Applications	ELISA, ID,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Diphtheria Toxin/Toxoid
----------	-------------------------