



Anti-DMPO polyclonal antibody (DPAB-DC4183)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	5,5-dimethyl-1-pyrroline-N-oxide (DMPO), octanoic acid adduct conjugated with ovalbumin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC, ELISA, EIA,
Format	Liquid
Size	100 µl
Buffer	In unpurified antiserum
Preservative	None
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DMPO
----------	------