



Anti-DNMT1 (N-terminal) polyclonal antibody (DPAB-DC4184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is specific to DNMT1.
Immunogen	N-terminus of DNMT1. The sequence is N-terminus
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), ICC, IF, IP,
Format	Liquid
Size	100 µl
Buffer	In antiserum
Preservative	None
Storage	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DNMT1
----------	-------

GENE INFORMATION

Entrez Gene ID [1786](#)

UniProt ID [I6L9H2](#)
