



Anti-DNP polyclonal antibody (DPAB-DC4185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This goat polyclonal antibody is specific to DNP and is intended for use in experiments related to receptor-mediated endocytosis and oxidative stress.
Immunogen	DNP conjugated with KLH.
Source/Host	Goat
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DNP
-----------------	-----