



Anti-Dopamine polyclonal antibody (DPAB-DC4188)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Dopamine conjugated with carrier protein.
Source/Host	Rat
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ICC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Dopamine (contracted from 3,4-dihydroxyphenethylamine) is a hormone and neurotransmitter of the catecholamine and phenethylamine families that plays a number of important roles in the human brain and body. Its name derives from its chemical structure: it
Keywords	Dopamine