



Anti-DCV Glycoprotein 1 Polyclonal antibody (DPAB-DC4190)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	DCV Glycoprotein 1
Immunogen	A synthetic peptide corresponding to C-terminus of Drosophila C virus DCVgp1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	DCV
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Drosophila C virus belongs to the genus Cripavirus and was previously thought to be a member of the virus family Picornaviridae; it has since been classified as belonging to the Dicistroviridae. It is a single stranded positive sense RNA virus of approxim
---------------------	---

Keywords

DCVgp1; Drosophila C virus DCVgp1
