



## Anti-DCV Glycoprotein 2 Polyclonal antibody (DPAB-DC4191)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	DCV Glycoprotein 2
<b>Immunogen</b>	A synthetic peptide corresponding to C-terminus of Drosophila C virus DCVgp2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	DCV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

### BACKGROUND

**Introduction** Drosophila C virus belongs to the genus Cripavirus and was previously thought to be a member of the virus family Picornaviridae; it has since been classified as belonging to the Dicistroviridae. It is a single stranded positive sense RNA virus of approxim

**Keywords**

DCVgp2; Drosophila C virus DCVgp2

---