



Anti-Dynein polyclonal antibody (DPAB-DC4197)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to rat Dynein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Dynein is a motor protein (also called molecular motor or motor molecule) in cells which converts the chemical energy contained in ATP into the mechanical energy of movement. Dynein transports various cellular cargo by "walking" along cytoskeletal microtu
Keywords	Dynein