



Anti-EGFP polyclonal antibody (DPAB-DC4198)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is known to react with EGFP confirmed by IHC and WB.
Immunogen	Recombinant protein corresponding to EGFP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	EGFP
----------	------