



Anti-EIA polyclonal antibody (DPAB-DC4199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant GST fusion protein corresponding to EIA.
Source/Host	Sheep
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB (Cell lysate), IP,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	EIA
-----------------	-----