



# Anti-procine neutrophil ELA2 (native) polyclonal antibody (DPAB-DC4200)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native purified ELA2 from procine neutrophil.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Pig
Purification	Affinity purification
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	ELA2
----------	------

## GENE INFORMATION

Entrez Gene ID [63036](#)

---

UniProt ID [P08217](#)

---