



Anti-Gamma-endorphin polyclonal antibody (DPAB-DC4204)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to endorphin-gamma. The sequence is CYGGFMTSEKSQTPLVTL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Gamma-endorphin is a polypeptide that is characterized by the presence of 17 amino acids, only one more than alpha-endorphin. Similar to other endorphins, research focusing upon gamma-endorphin has been ongoing since its discovery in the 1970s. Yet, most
Keywords	Gamma-endorphin; Endorphin-gamma