



Anti-*E. cloacae* Polyclonal antibody (DPAB-DC4205)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>E. cloacae</i>
Immunogen	<i>Enterobacter cloacae</i> total protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>E. cloacae</i>
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction *Enterobacter cloacae* is a clinically significant Gram-negative, facultatively-anaerobic, rod-

shaped bacterium.

Keywords

E. cloacae; Enterobacter cloacae
