



Anti-E. faecium Polyclonal antibody (DPAB-DC4206)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Specificity | Whole cells. This antiserum has not been absorbed and may react with related microorganisms. |
| Target | E. faecium |
| Immunogen | Enterococcus faecium |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | E. faecium |
| Conjugate | Unconjugated |
| Applications | IF, |
| Format | Liquid |
| Size | 100 µl |
| Buffer | In 10 mM PBS, pH 7.2 (0.1% sodium azide) |
| Preservative | 0.1% Sodium Azide |
| Storage | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

Introduction

Enterococcus faecium is a Gram-positive, alpha hemolytic or nonhemolytic bacterium in the genus *Enterococcus*. It can be commensal (innocuous, coexisting organism) in the human intestine, but it may also be pathogenic, causing diseases such as neonatal men

Keywords

E. faecium; *Enterococcus* species